



# Piloter Vet

A Smart Tablet Color Doppler Ultrasound System

**wisonic**

# Vet Application

Pet Hospitals and Clinics

1. General application
2. Sterilization

Horse Riding Club

1. Physical examination
2. Breeding and monitoring equine pregnancy
3. Tendon scanning

Livestock

1. Bovine, equine, ovine pregnancy
2. Musculoskeletal
3. Meat and fat quality rating

Aquaculture

1. Gender screening
2. Breeding selection

Safari & Zoo

1. Regular physical examination
2. Reproduction



Table Stand

## Dedicated Transducers



Convex C5-2  
Medium to large animals  
Abdomen, reproductive system



Linear L15-4NB  
MSK of medium to large animals  
Small animals whole body



Micro-Convex C8-3  
Small animals abdomen and heart




Phased P7-3  
Small to medium animals heart




Linear L9-4V  
Intracavity, reproductive system


# Excellent Image Quality




Liver




Splenic Nodule




Bladder Wall



Intestine




Dog Cardiac



Kidney Blood Flow


# Broad-Range of Transducers

Piloter is capable of performing versatile applications from small animals to large animals with excellent image quality. The ergonomic button probe can adjust gain, depth, and freeze/unfreeze the image, makes diagnosis easier.



# wiCloud

Remote diagnosis and remote control assist veterinaries to diagnose more confidently and efficiently



# QR Code

Help veterinaries to share and communicate cases easily



## Shenzhen Wisonic Medical Technology Co.,Ltd.

Address: 5th Floor, No.6 Building, Pingshan Tech Park,  
Taoyuan Street , Nanshan, Shenzhen, 518055, P.R.China.

Tel: +86-755-86007788

Fax: +86-755-86007799

Email: [info@wisonic.com](mailto:info@wisonic.com)

[www.wisonic.com](http://www.wisonic.com)

© Shenzhen Wisonic Medical Technology Co., Ltd. All rights reserved. **wisonic** is a registered trademark of Shenzhen Wisonic Medical Technology Co., Ltd. The internal information is for learning only. Shenzhen Wisonic Medical Co., Ltd. reserves the right to modify product specifications and / or terminate the production of any product at any time. Product specifications are subject to change without notice.