




D-Pro
Veterinary Dental Unit

RWD



Every veterinary dentist
deserves a **professional
animal dental unit** 

Buy one dental unit, you will get ...

- ✓ Shortening the treatment and anesthesia to reduce the risk of anesthesia.
- ✓ Reducing the damage of alveolar bone and oral mucosa caused by violent tooth extraction.
- ✓ Reducing the risk of retained roots.
- ✓ More flexible dental treatment programs.
- ✓ High return on equipment investment.

D-pro veterinary dental unit provides a variety of functions for veterinarians, such as grinding, polishing, cutting and cleaning, and makes daily animal oral health work easy and efficient.



Water and electricity separation design to ensure the safety of the equipment.



360 ° hand piece, light in weight, can reduce hand fatigue effectively.



400000 rpm high-speed LED hand piece, weighing only 66g, can help veterinarians easily cut animal teeth.



5000 rpm low-speed hand piece helps the veterinarian to further clean the plaque and inhibit the growth of the plaque.



Oil free silent air compressor with noise of 60dB, provides a comfortable working environment for veterinarians.



With a 20L large capacity air tank to ensure the lasting supply.



There is a magnetic platform to absorb the burs and other fine dental tools to avoid accidental sliding or loss.



The whole machine equipped with universal wheel to ensure flexible movement.

Configuration

Item	Specification	
Out put pressure of air compressor	85psi (6kgf/cm ²)	
High speed hand piece	Working pressure	40psi (2.8kgf/cm ²)
	Base speed	400000±10%rpm
	Weight	66g
Low speed hand piece	Working pressure	35-40psi
	Base speed	5000±10%rpm
	Weight	95g
Air tank	20L	

RWD

RWD Life Science Co.,Ltd

Add: 2F, Robeta Building, No.11 High Tech North 4th Road, Nanshan District, Shenzhen City, Guangdong Province, 518057, China.
Tel: (86)-755-86111286 Web: www.rwdls.com E-mail: rwd@rwdmall.com

RWD Life Science Inc.

Add: 6540 Lusk Blvd, San Diego, CA92121, USA
Tel: 001-858-9006602 Web: www.rwdstco.com E-mail: market@rwdstco.com