

Vet HF

Digital Electro Surgery Units

Excellent performance • Excellent price

Vet HF-M

Microprocessor-controlled monopolar electro surgery unit

Safe and reliable device

For little or medium surgery



Features:

- Maximal output power: 100 Watt
- 4 operation modes:
 1. CUT - smooth cutting mode
 2. CUT1 - cutting with heightened degree of homeostasis
 3. COAG1 - forced (standard) coagulation
 4. COAG2 - soft coagulation
- Independent control cut / coagulation
- Digital display from actual setting
- Automatic saving of the last settings
- Control with a foot switch
- Sensor for alarm indication in case of an overheating or a technical defective, automatic control of the neutral electrode
- Auto-Diagnostic
- Acoustic and visual alarm

Technical Data:

| | |
|--------------------------|--|
| Accessories: | Footswitch, reusable handle, 7 active electrodes, neutral plate and main power cable |
| Working frequency: | 500 KHz |
| Modulation frequency: | 33 KHz |
| Activation: | Footswitch |
| Dimensions (L x W x H) : | 255 x 215 x 85 mm |
| Weight: | 5 kg |


Vet HF surgery units

Vet HF-B

Microprocessor-controlled bipolar electrosurgery unit

Safe and reliable device

Perfect for the small animal practice



Features:

- Maximal output power: 100 Watt
- Operation modes:
 1. CUT - smooth cutting mode
 2. CUT1 - cutting with heightened degree of haemostasis
 3. COAG1 - forced (standard) coagulation
 4. COAG2 - soft coagulation
 5. Bipolar coagulation
 6. Bipolar micro-coagulation
- Independent control cut / coagulation
- Digital display from actual setting
- Automatic saving of the last settings
- Control with a foot or a hand switch
- Sensor for alarm indication in case of an overheating or a technical defective, automatic control of the neutral electrode
- Auto-Diagnostic
- Acoustic and visual alarm

Technical Data:

| | |
|--------------------------|--|
| Accessories: | Twin footswitch with yellow and blue sections, straight bipolar forceps, 7 active electrodes, neutral plate and main power cable |
| Working frequency: | 500 KHz |
| Modulation frequency: | 33 KHz |
| Activation: | Twin footswitch and Hand Switched Pencil |
| Dimensions (L x W x H) : | 255 x 215 x 85 mm |
| Weight: | 5 kg |



Medical equipment for veterinarians

15, rue du Noroit
22100 Taden - France

Tél. +33 (0)2 96 85 86 76

Fax +33 (0)2 96 85 86 77

welcome@manomedical.com

www.manomedical.com