

7. TECHNICAL DATA SHEET

Product Name	Stabilo Compact 50x100 Stabilo Compact 60x100	Stabilo Super with or without wheels	Stabilo Original with or without wheels	Stabilo Elite with or without wheels
Product Code	STABI02 / STABI05	STABI0 / STABI1	STABI03 / STABI04	STABI06 / STABI07
Safe working load	200kg	200kg	200kg	250kg
Max Animal weight	175kg	175kg	175kg	225kg
Height adjustment:	45 to 100cm with a speed of 7,5mm/s - 8.000N, 470W, 2.0Ah, stroke 300mm	38 to 97cm with a speed of 7,5mm/s - 8.000N, 470W, 2.0Ah, stroke 300mm	38 to 97cm with a speed of 7,5mm/s - 8.000N, 470W, 2.0Ah, stroke 200mm	38 to 110cm with a speed of 7,5mm/s - 8.000N, 470W, 2.0Ah, stroke 300mm
Adjustment system	German electric actuator	German electric actuator	German electric actuator	German electric actuator
Adjusting electrics	Foot- or handswitch	Foot- or handswitch	Foot- or handswitch	Foot- or handswitch
Wheels	Ø75mm	Ø75mm	Ø75mm	Ø75mm
Braking system	Centralized control on the back of the table	Centralized control on the back of the table	Centralized control on the back of the table	Centralized control on the back of the table
Equipment weight	± 55kg	± 70kg	± 75kg	± 80kg
Length	100cm	120cm	120cm	120cm
Width	50cm / 60cm	65cm	65cm	65cm
Workplate thickness	20mm	20mm	20mm	20mm
Filling	Wood	Wood	Wood	Wood
Top	Grey diamond plate rubber	Grey diamond plate rubber	Grey diamond plate rubber	Grey diamond plate rubber